
pybdf Documentation

Release alpha3

Samuele Carcagno

July 16, 2012

CONTENTS

1	Introduction	3
2	Download and Installation	5
2.1	Download	5
2.2	Installation	5
3	Usage	7
4	Bugs	9
5	Benchmarks	11
6	pybdf – Class to read BIOSEMI BDF files	13
7	Indices and tables	17
	Python Module Index	19
	Index	21

Contents:

INTRODUCTION

Author Samuele Carcagno

pybdf is a pure python library to read BIOSEMI 24-bit BDF files. While being slower than alternative C-based libraries like *BioSig* <<http://biosig.sourceforge.net/>>, it is very easy to install and use.

DOWNLOAD AND INSTALLATION

2.1 Download

The latest release of pybdf can be downloaded from the python package index:

<http://pypi.python.org/pypi/pybdf/>

For developers: the source code of pybdf is hosted on launchpad:

<https://launchpad.net/pybdf>

2.2 Installation

2.2.1 Requirements

To install pybdf you will need:

- python \geq 3.2
- numpy

On all platforms, after having unpacked the archive you can install pybdf by running:

```
python setup.py install
```

On Windows, you can alternatively use the binary installer, if provided. Note that pybdf has not been extensively tested on Windows.

USAGE

To open a bdf file you need to create a `bdfRecording` object as follows:

```
bdf_rec = bdfRecording('res1.bdf')
```

you can then query the properties of therecording stored in the BDF header using the appropriate functions, which are fully described [here](#). Some examples are shown below.

Get the duration of the recording:

```
bdf_rec.recordDuration
```

Get the sampling rate of each channel:

```
bdf_rec.sampRate
```

Get the channel labels:

```
bdf_rec.chanLabels
```

There are two functions to read in the data. The first function reads each channel sequentially:

```
rec = bdf_rec.get_data()
```

the second function reads the channels in parallel, and is thus faster on multicore machines:

```
rec = bdf_rec.get_data_parallel()
```

either function returns the same result, that is a python dictionary with the following fields: - `data` : an array of floats with dimenions `nChannels X nDataPoints` - `trigChan` : an array of integers with the triggers in decimal format - `statusChan` : an array of integers with the status codes in decimal format For example, to get the value of the first sample of the recording, in the first channel, you can type:

```
rec['data'][0,0]
```

the same sample value, but for the second channel, is stored in:

```
rec['data'][1,0]
```

`trigChan` contains the triggers for the experimental conditions, in decimal format. The `statusChan`, on the other hand, contains system codes, like `cm in/out-of range`, `battery low/OK`.

Other usage examples are provided in the 'examples' directory inside the `pybdf` source archive.

Beware that `pybdf` does not check that you have sufficeint RAM to read all the data in a bdf file. If you try to read a file that is too big for your hardware, you system may become slow or unresponsive. Initially try reading only a small amount of data, and check how much RAM that uses. You can read only a portion of the data by passing the beginning

and end arguments to the `get_data()` or `get_data_parallel()` functions. For example, to read the first 10 seconds of the recording, use:

```
rec = bdf_rec.get_data_parallel(beginning=0, end=10)
```

BUGS

Please, report any bugs on Launchpad <https://launchpad.net/pybdf>

Currently there are problems with the `get_data_parallel()` function on Windows. Please, use the `get_data()` function instead, or use Linux.

BENCHMARKS

To give you an idea of the speed of pybdf, here are some rough benchmarks.

Using the `get_data_parallel()` function:

Channels	Duration (s)	Samp. Rate	File (MB)	CPU	Time (s)
9	931	2048	49.1	Intel Core i7-870	10
9	1457	2048	76.8	Intel Core i7-870	15
41	651	2048	156.4	Intel Core i7-870	21
9	931	2048	49.1	Intel Core2 Quad Q6600	16
9	1457	2048	76.8	Intel Core2 Quad Q6600	24
41	651	2048	156.4	Intel Core2 Quad Q6600	31

PYBDF – CLASS TO READ BIOSEMI BDF FILES

This module can be used to read the header and data from 24-bit BIOSEMI BDF files recorded with the ActiveTwo system.

```
>>> bdf_rec = bdfRecording('res1.bdf') #create bdfRecording object
>>> bdf_rec.recordDuration #how many seconds the recording lasts
>>> bdf_rec.sampRate #sampling rate for each channel
>>> #read 10 seconds of data from the first two channels
>>> rec = bdf_rec.get_data(channels=[0, 1], beginning=0, end=10)
>>> rec = bdf_rec.get_data_parallel() #read all data using multiprocessing
```

class `pybdf.bdfRecording` (*fileName*)

Class for dealing with BIOSEMI 24-bit BDF files. A `bdfRecording` object is created with the following syntax:

```
>>> bdf_rec = bdfRecording('bdf_file.bdf')
```

This reads the BDF header, but not the data. You need to use the `get_data()` or `get_data_parallel()` methods to read the data. The full documentation of the BDF file format can be found here: http://www.biosemi.com/faq/file_format.htm

idCode [str] Identification code

subjId [str] Local subject identification

recId [str] Local recording identification

startDate [str] Recording start date

startTime [str] Recording start time

nBytes [int] Number of bytes occupied by the bdf header

versionDataFormat [str] Version of data format

nDataRecords [int] Number of data records “-1” if unknown

recordDuration [float] Duration of a data record, in seconds

nChannels [int] Number of channels in data record

chanLabels [list of str] Channel labels

transducer [list of str] Transducer type

physDim [str] Physical dimension of channels

physMin [list of int] Physical minimum in units of physical dimension

physMax [list of int] Physical maximum in units of physical dimension

digMin [list of int] Digital minimum

digMax [list of int] Digital maximum

prefilt [list of str] Prefiltering

nSampRec [list of int] Number of samples in each data record

reserved [list of str] Reserved

scaleFactor [list of floats] Scaling factor for digital to physical dimension

sampRate [list of int] Recording sampling rate

statusChanIdx [int] Index of the status channel

nDataChannels [int] Number of data channels containing data (rather than trigger codes)

dataChanLabels [list of str] Labels of the channels containing data (rather than trigger codes)

get_data (*beginning=0, end=None, channels=None, trig=True, status=True, norm_trig=True, norm_status=True*)
Read the data from a bdfRecording object

beginning [int] Start time of data chunk to read (seconds).

end [int] End time of data chunk to read (seconds).

channels [list of integers or strings] Channels to read. Both channel numbers, or channel names are accepted. Note that channel numbers are indexed starting from *zero*.

trig [boolean] If True, return the channel containing the triggers

status [boolean] If True, return the channel containing the status codes

norm_trig [boolean] If True, the trigger channel will only signal *changes* between one trigger status to the next. A trigger value that is equal to the previous one will be set to zero

norm_status [boolean] If True, the status channel will only signal *changes* between one status code to the next. A code value that is equal to the previous one will be set to zero

rec [a dictionary with three keys]

- **data** : an array of floats with dimensions nChannels X nDataPoints
- **trigChan** : an array of integers with the triggers in decimal format
- **statusChan** : an array of integers with the status codes in decimal format
- **chanLabels** : a list containing the labels of the channels that were read, in the same order they are inserted in the data matrix

```
>>> x = bdfRecording('res1.bdf')
>>> rec = x.get_data(channels=[0, 2], beginning=0, end=10)
```

get_data_parallel (*beginning=0, end=None, channels=None, trig=True, status=True, norm_trig=True, norm_status=True*)
Read the data from a bdfRecording object using the multiprocessing module to exploit multicore machines.

beginning [int] Start time of data chunk to read (seconds).

end [int] End time of data chunk to read (seconds).

channels [list of integers or strings] Channels to read. Both channel numbers, or channel names are accepted. Note that channel numbers are indexed starting from *zero*.

trig [boolean] If True, return the channel containing the triggers

status [boolean] If True, return the channel containing the status codes

norm_trig [boolean] If True, the trigger channel will only signal *changes* between one trigger status to the next. A trigger value that is equal to the previous one will be set to zero

norm_status [boolean] If True, the status channel will only signal *changes* between one status code to the next. A code value that is equal to the previous one will be set to zero

rec [a dictionary with three keys]

- **data** : an array of floats with dimensions nChannels X nDataPoints
- **trigChan** : an array of integers with the triggers in decimal format
- **statusChan** : an array of integers with the status codes in decimal format
- **chanLabels** : a list containing the labels of the channels that were read, in the same order they are inserted in the data matrix

```
>>> x = bdfRecording('res1.bdf')
>>> rec = x.get_data_parallel(channels=[0, 2], beginning=0, end=10)
```


INDICES AND TABLES

- *genindex*
- *modindex*
- *search*

PYTHON MODULE INDEX

p

pybdf, 13

INDEX

B

bdfRecording (class in pybdf), 13

G

get_data() (pybdf.bdfRecording method), 14

get_data_parallel() (pybdf.bdfRecording method), 14

P

pybdf (module), 13